



**PM200**

**PM100A/D/USB**

**PM160**

**Write Your Own Application**

- **Drivers**
  - **Samples**
  - **SCPI Command Reference**
-

# Contents

<b>Foreword</b>	<b>0</b>
<b>1 Write Your Own Application</b>	<b>2</b>
1.1 Windows XP 32bit _____	3
1.2 Windows Vista / 7 / 8 - 32bit _____	4
1.3 Windows Vista / 7 / 8 - 64bit _____	5
1.4 Simple LabVIEW Example using SCPI commands _____	7
<b>2 SCPI Commands</b>	<b>11</b>
2.1 An Introduction to the SCPI language _____	11
2.2 IEEE488.2 Common Commands _____	15
2.2.1 Command Summary _____	15
2.2.2 Command Reference _____	16
2.2.3 PM200 specific SCPI Command Reference _____	17
2.2.3.1 SYSTem subsystem commands _____	17
2.2.3.2 STATus subsystem commands _____	18
2.2.3.3 DISPlay subsystem commands _____	19
2.2.3.4 CALibration subsystem commands _____	19
2.2.3.5 SENSE subsystem commands _____	19
2.2.3.6 INPut subsystem commands _____	22
2.2.3.7 Measurement commands _____	23

---

# 1 Write Your Own Application

In order to write your own application, you need a specific instrument driver and some tools for use in different programming environments. The driver and tools are included in the installer package and cannot be found as separate files on the installation CD.

In this section the location of drivers and files, required for programming in different environments, are given for installation under Windows XP (32 bit) and Windows 7 (32 and 64 bit)

## Note

PM200 software and drivers contains 32 bit and 64 bit applications.

In 32 bit systems, only the 32 bit components are installed to

**C:\Program Files\...**

In 64 bit systems the 64 bit components are being installed to

**C:\Program Files\...**

while 32 bit components can be found at

**C:\Program Files (x86)\...**

In the table below you will find a summary of what files you need for particular programming environments.

Programming environment	Necessary files
<b>C, C++, CVI</b>	*.fp (function panel file; CVI IDE only) *.h (header file) *.lib (static library) *.dll (dynamic linked library)
<b>C#</b>	.net wrapper dll
<b>Visual Studio</b>	*.h (header file) *.lib (static library) or .net wrapper dll
<b>LabView</b>	*.fp (function panel) and NI VISA instrument driver Beside that, LabVIEW driver vi's are provided with the *.llb container file

## Note

All above environments require also the NI VISA instrument driver dll !

During NI-VISA Runtime installation, a system environment variable `VXIPNPPATH` for including files is created. It contains the information where the drivers are installed to, usually to `C:\Program Files\IVI Foundation\VISA\WinNT\`.

This is the reason, why after installation of a NI-VISA Runtime a system reboot is required: This environment variable is necessary for installation of the instrument driver software components.

In the next sections the location of above files is described in detail.

## 1.1 Windows XP 32bit

### NI VISA Instrument Driver 32bit

C:\Program Files\IVI Foundation\VISA\WinNT\Bin\PM100D\_32.dll

### Online Help for NI VISA Instrument Driver

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\Manual\PM100D.html

### NI LabVIEW Driver

The LabVIEW version of this driver was generated with the "NI LabVIEW Instrument Driver Import Wizard 2.0" in conjunction with LabVIEW 2011.

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\LabVIEW\PM100D.llb

### Header File

C:\Program Files\IVI Foundation\VISA\WinNT\include\PM100D.h

### Static Library

C:\Program Files\IVI Foundation\VISA\WinNT\lib\msc\PM100D\_32.lib

### Function Panel

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\PM100D.fp

### .net wrapper dll

C:\Program Files\IVI Foundation\VISA\VisaCom\...  
...Primary Interop Assemblies\Thorlabs.PM100D.dll

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\DotNet\...  
...Thorlabs.PM100D.dll

C:\Program Files\Microsoft.NET\Primary Interop Assemblies\...  
...Thorlabs.PM100D.dll

## Examples

### ANSI-C

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\Samples\C\sample.c

Thorlabs PM100x/PM160/PM200 Instrument Driver Sample Application with console interface. Read the comment text in the sample.c file for more information.

### C# Visual Studio 2010

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\Samples\...  
...DotNet\DotNetSample.csproj

Very basic C# sample project. The application opens a session to a device takes one measurement and closes the device again.

## 1.2 Windows Vista / 7 / 8 - 32bit

### NI VISA Instrument Driver

C:\Program Files\IVI Foundation\VISA\WinNT\Bin\PM100D\_32.dll

### Online Help for NI VISA Instrument Driver:

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\Manual\PM100D.html

### NI LabVIEW Driver

The LabVIEW version of this driver was generated with the "NI LabVIEW Instrument Driver Import Wizard 2.0" in conjunction with LabVIEW 2011.

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\LabVIEW\PM100D.llb

### Header File

C:\Program Files\IVI Foundation\VISA\WinNT\include\PM100D.h

### Static Library

C:\Program Files\IVI Foundation\VISA\WinNT\lib\msc\PM100D\_32.lib

### Function Panel

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\PM100D.fp

### .net wrapper dll

C:\Program Files\IVI Foundation\VISA\VisaCom\...

...Primary Interop Assemblies\Thorlabs.PM100D.dll

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\DotNet\...

...Thorlabs.PM100D.dll

C:\Program Files\Microsoft.NET\Primary Interop Assemblies\...

...Thorlabs.PM100D.dll

### Examples:

#### ANSI-C

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\Samples\C\sample.c

Thorlabs PM100x/PM160/PM200 Driver Sample Application with console interface. Read the comment text in the sample.c file for more information.

#### C# Visual Studio 2010

C:\Program Files\IVI Foundation\VISA\WinNT\PM100D\Samples\...

...DotNet\DotNetSample.csproj

Very basic C# sample project. The application opens a session to a device takes one measurement and closes the device again.

## 1.3 Windows Vista / 7 / 8 - 64bit

### NI VISA Instrument Driver 32bit

C:\Program Files (x86)\IVI Foundation\VISA\WinNT\Bin\PM100D\_32.dll

### NI VISA Instrument Driver 64bit

C:\Program Files\IVI Foundation\VISA\Win64\Bin\PM100D\_64.dll

### Online Help for NI VISA Instrument Driver

C:\Program Files\IVI Foundation\VISA\Win64\PM100D\Manual\PM100D.html

### NI LabVIEW Driver 32bit

C:\Program Files (x86)\IVI Foundation\VISA\WinNT\PM100D\...  
...LabVIEW\PM100D.llb

### NI LabVIEW Driver 64bit

C:\Program Files\IVI Foundation\VISA\Win64\PM100D\LabVIEW\PM100D.llb

### Header File 32bit

C:\Program Files (x86)\IVI Foundation\VISA\WinNT\include\PM100D.h

### Header File 64bit

C:\Program Files\IVI Foundation\VISA\Win64\include\PM100D.h

### Static Library 32bit

C:\Program Files (x86)\IVI Foundation\VISA\WinNT\lib\msc\PM100D\_32.lib

### Static Library 64bit

C:\Program Files\IVI Foundation\VISA\Win64\Lib\_x64\msc\PM100D\_64.lib

### Function Panel 32bit

C:\Program Files (x86)\IVI Foundation\VISA\WinNT\PM100D\PM100D.fp

### Function Panel 64bit

C:\Program Files\IVI Foundation\VISA\Win64\PM100D\PM100D.fp

**.net wrapper dll 32bit**

```
C:\Program Files (x86)\IVI Foundation\VISA\VisaCom\...  
...Primary Interop Assemblies\Thorlabs.PM100D.dll  
C:\Program Files (x86)\Microsoft.NET\Primary Interop Assemblies\...  
...Thorlabs.PM100D.dll
```

**.net wrapper dll 64bit**

```
C:\Program Files\IVI Foundation\VISA\VisaCom64\...  
...Primary Interop Assemblies\Thorlabs.PM100D.dll
```

**Examples:****ANSI-C**

```
C:\Program Files\IVI Foundation\VISA\Win64\PM100D\Samples\C\sample.c
```

Thorlabs PM100x/PM160/PM200 Driver Sample Application with console interface. Read the comment text in the sample.c file for more information.

**C# Visual Studio 2010**

```
C:\Program Files\IVI Foundation\VISA\Win64\PM100D\Samples\DotNet\...  
...DotNetSample_64.csproj
```

Very basic C# sample project. The application opens a session to a device takes one measurement and closes the device again.

## 1.4 Simple LabVIEW Example using SCPI commands


The PM200 Instrument driver is not required for this LabVIEW example.

PM100D  
Simple  
Example

### PM100D Simple Example.vi


This VI shows how to communicate with a PM200 optical power/energy meter with SCPI commands. The following steps are demonstrated within this application:

- Initializing the instrument
- Getting system info
- Setting parameters
- Measurement configuration
- Measuring queue
- Fetching and displaying a measurement value
- Closing the application





## Block diagram



**Averaging Rate**

Sets the averaging rate - 3000 averages take approximately 1s for performing a new measurement value

**Readout Config**

Selects the measurement parameter

- POW power measurement in W
- CURR current measurement in A
- VOLT voltage measurement in V
- ENER energy measurement in J
- FREQ frequency / repetition rate measurement in Hz
- PDEN power density measurement in W/cm<sup>2</sup>
- EDEN energy density measurement in J/cm<sup>2</sup>
- RES thermistor resistance measurement in Ohm
- TEMP temperature measurement in °C

**Stop**

Stops application

**Timeout [ms]**

Sets a timeout value in ms that allows the instrument to sample.

The timeout must be longer than it takes to perform a new measurement. This has especially to be considered when performing single shot energy measurements.

**Error log**

Error indicator, 'no error' is suppressed

**Readout string**

Formatted measurement value

- limitation to 4 significant digits
- SI notation
- '.' decimal separator
- appended unit according readout configuration

**New**

Indicator lights up, when a new measurement value is processed

**Sensor Flags**

Sensor flag bitmap:

- 1 Is power sensor
- 2 Is energy sensor
- 16 Response settable
- 32 Wavelength settable
- 64 Tau settable
- 256 Has temperature sensor

**Sensor Name**

Name of connected power/energy sensor

**Sensor SN**

Serial number of connected power/energy sensor

**IDN String**

Answer from instrument on \*idn? command:

- manufacturer

- model number
- serial number
- firmware version

**Readout value**

Plain readout value in full resolution

**PM100D\_Initialize.vi**

This VI scans for connected devices, that can be selected in a dialog box. Next steps

- setting timeout
- performing identification query
- configuring the operation register to '512'; flag gets to HI when a new measurement value is ready to fetch
- clear operation register

**PM100D\_Write.vi**

Writes a SCPI command to the connected instrument

**PM100D\_Read.vi**

Reads data from the connected instrument. All query commands according the SCPI command table are terminated by a question mark (?)

**PM100D\_SensorFlag.vi**

This VI queries all relevant sensor info with SYST:SENS:IDN?

- sensor name
- sensor serial number
- calibration message
- sensor type
- sensor sub-type
- sensor flags

**PM100D\_ConfMeas.vi**

- CONFigure measurement; CONF:<parameter> (POW, ENER, etc.)
- ABORt measurement
- Clear operation register with STAT:OPER?
- INITiate measurement

**PM100D\_SYST-ERR.vi**

This VI lists all errors coming from the instrument with the command SYST:ERR?  
'no error' is suppressed

**PM100D\_Close.vi**

Closes the VISA session

Sets the connected instrument in local mode (default option)

## 2 SCPI Commands

This section describes in detail the SCPI command set with respect to PM200 power meter console. However, these commands are compatible with PM100A, PM100D and PM100USB consoles as well as with the PM160 Hand-held Power Meter, except some console and sensor related functions.

### 2.1 An Introduction to the SCPI language

The PM200 interface commands use the SCPI (Standard Commands for Programmable Instruments), an ASCII-based command language that was designed for test and measurement instruments.

SCPI commands are based on a hierarchical structure, also known as a tree system. In this system, associated commands are grouped together under a common node or root, thus forming subsystems. A portion of the SENSE subsystem is shown below to illustrate the tree system.

#### SENSe:

```
CORRection
  :COLlect
    :ZERO
      [:INITiate]
      :ABORt
      :STATe?
      :MAGNitude?
:BEAMdiameter {MINimum|MAXimum|DEFault|<numeric_value>[mm]}
:BEAMdiameter? [{MINimum|MAXimum|DEFault}]
:WAVelength {MINimum|MAXimum|<numeric_value>[nm]}
:WAVelength? [{MINimum|MAXimum}]
:POWer
  [:PDIode]
    [:RESPonse] MINimum|MAXimum|DEFault|<numeric_value>[A]}
    [:RESPonse]? [{MINimum|MAXimum|DEFault}]
  :THERmopile
    [:RESPonse] {MINimum|MAXimum|DEFault|<numeric_value>[V]}
    [:RESPonse]? [{MINimum|MAXimum|DEFault}]
```

SENSe is the root keyword of the command, CORRection is the second-level keyword, and COLlect and BEAMdiameter are third-level keywords, and so on.

A colon (:) separates a command keyword from a lower-level keyword.

#### Command Format

The format used to show commands in this manual is shown below:

```
CURRent[:DC]:RANGe {MINimum|MAXimum|<numeric_value>[A]}
CORRection:BEAMdiameter {MINimum|MAXimum|DEFault|<numeric_value>[mm]}
```

The command syntax shows most commands (and some parameters) as a mixture of upper- and lower-case letters. The upper-case letters indicate the abbreviated spelling for the command. For shorter program lines, send the abbreviated form. For better program readability, send the long form.

For example, in the above syntax statement, `CURR` and `CURRENT` are both acceptable forms. You can use upper- or lower-case letters. Therefore, `CURRENT`, `current` and `Current` are all acceptable. Other forms, such as `CUR=` and `CURREN`, will generate an error.

**Braces ( { } )** enclose the parameter choices for a given command string. The braces are not sent with the command string. A *vertical bar* ( | ) separates multiple parameter choices for a given command string.

**Triangle brackets ( < > )** indicate that you must specify a value for the enclosed parameter. For example, the above syntax statement shows the *range* parameter enclosed in triangle brackets. The brackets are not sent with the command string. You must specify a value for the parameter (such as `"CURR:DC:RANG 50E-6"`).

Some parameters are enclosed in *square brackets* ( [ ] ). The brackets indicate that the parameter is optional and can be omitted. The brackets are not sent with the command string. In this example `[:DC]` can be omitted, so the command string can be shortened to `"CURR:RANG 50E-6"`. If you do not specify a value for an optional parameter, the power/energy meter chooses a default value.

### Command Separators

A *colon* ( : ) is used to separate a command keyword from a lower-level keyword. You must insert a *blank space* to separate a parameter from a command keyword. If a command requires more than one parameter, you must separate adjacent parameters using a *comma* as shown below:

```
"SYST:TIME 10, 34, 48"
```

A *semicolon* ( ; ) is used to separate commands within the *same* subsystem, and can also minimize typing. For example, sending the following command string:

```
"CORR:BEAM 1; WAV 1310"
```

... is the same as sending the following two commands:

```
"CORR:BEAM 1"  
"CORR:WAV 1310"
```

Use a colon and a semicolon to link commands from different subsystems. For example, in the following command string, an error is generated if you do not use both the colon and semicolon:

```
"CORR:BEAM 1;:AVER 300"
```

### Using the *MIN* and *MAX* Parameters

You can substitute `MINimum=` or `MAXimum=` in place of a parameter for many commands. For example, consider the following command:

```
CURRent[:DC]:RANGe {MINimum|MAXimum|<numeric_value>[A]}
```

Instead of selecting a specific current range, you can substitute `MIN` to set the range to its minimum value or `MAX` to set the range to its maximum value.

### Querying Parameter Settings

You can query the current value of most parameters by adding a question mark ( ? ) to the command. For example, the following command sets the operating wavelength to 1550 nm:

```
"CORR:WAV 1550"
```

You can query the operating wavelength by executing: `"CORR:WAV?"`

You can also query the minimum or maximum operating wavelength allowed as follows:

```
"CORR:WAV? MIN"
```

```
"CORR:WAV? MAX"
```

**Caution**

*If you send two query commands without reading the response from the first, and then attempt to read the second response, you may receive some data from the first response followed by the complete second response. To avoid this, do not send a query command without reading the response. When you cannot avoid this situation, send a device clear before sending the second query command.*

**SCPI Command Terminators**

A command string sent to the power/energy meter must terminate with a <new line> character. The IEEE-488 EOI (end-or-identify) message is interpreted as a <new line> character and can be used to terminate a command string in place of a <new line> character. A <carriage return> followed by a <new line> is also accepted. Command string termination will always reset the current SCPI command path to the root level.

**IEEE488.2 Common Commands**

The IEEE-488.2 standard defines a set of common commands that perform functions like reset, self-test, and status operations. Common commands always begin with an asterisk (\*), are four to five characters in length, and may include one or more parameters. The command keyword is separated from the first parameter by a blank space. Use a semicolon (;) to separate multiple commands as shown below:

```
"*RST; *CLS; *ESE 32; *OPC?"
```

**SCPI Parameter Types**

The SCPI language defines several different data formats to be used in program messages and response messages.

**Numeric Parameters** Commands that require numeric parameters will accept all commonly used decimal representations of numbers including optional signs, decimal points, and scientific notation.

Special values for numeric parameters like `MINimum`, `MAXimum` and `DEFault` are also accepted. You can also send engineering unit suffixes with numeric parameters (e.g., M, K, or u). If only specific numeric values are accepted, the power/energy meter will automatically round the input numeric parameters. The following command uses a numeric parameter:

```
POWer:REFeRence {MINimum|MAXimum|DEFault|<numeric_value>[W]}
```

**Discrete Parameters** Discrete parameters are used to program settings that have a limited number of values (like W, DBM). They can have a short form and a long form just like command keywords. You can mix upper- and lower-case letters. Query responses will *always* return the short form in all upper-case letters. The following command uses discrete parameters:

```
POW:UNIT {W|DBM}
```

**Boolean Parameters** Boolean parameters represent a single binary condition that is either true or false. For a false condition, the power/energy meter will accept "OFF" or "0". For a true condition, the meter will accept "ON" or "1". When you query a boolean setting, the instrument will *always* return "0" or "1". The following command uses a boolean parameter:

```
CURRent:RANGe:AUTO {OFF|0|ON|1}
```

**String Parameters** String parameters can contain virtually any set of ASCII characters. A string *must* begin and end with matching quotes; either with a single quote or with a double quote. You can include the quote delimiter as part of the string by typing it twice without any characters in between. The following command uses a string parameter:

```
DIAG:CALString <quoted string>
```

## 2.2 IEEE488.2 Common Commands

Common commands are device commands that are common to all devices according to the IEEE488.2 standard. These commands are designed and defined by this standard. Most of the commands are described in detail in this section. The following common commands associated with the status structure are covered in the “Status Structure” section: \*CLS, \*ESE, \*ESE?, \*ESR?, \*SRE, \*SRE?, \*STB?

### 2.2.1 Command Summary

Mnemonic	Name	Description
<b>*CLS</b>	Clear status	Clears all event registers and Error Queue
<b>*ESE</b> <b>&lt;NRf&gt;</b>	Event enable command	Program the Standard Event Enable Register
<b>*ESE?</b>	Event enable query	Read the Standard Event Enable Register
<b>*ESR?</b>	Event status register query	Read and clear the Standard Event Register
<b>*IDN?</b>	Identification query	Read the unit's identification string
<b>*OPC</b>	Operation complete command	Set the Operation Complete bit in the Standard Event Register
<b>*OPC?</b>	Operation complete query	Places a “1” into the output queue when all device operations have been completed
<b>*RST</b>	Reset command	Returns the unit to the *RST default condition
<b>*SRE</b> <b>&lt;NRf&gt;</b>	Service request enable command	Programs the Service Request Enable Register
<b>*SRE?</b>	Service request enable query	Reads the Service Request Enable Register
<b>*STB?</b>	Status byte query	Reads the Status Byte Register
<b>*TST?</b>	Self-test query	Performs the unit's self-test and returns the result.
<b>*WAI</b>	Wait-to-continue command	Wait until all previous commands are executed



## 2.2.2 Command Reference

### **\*IDN? – identification query - read identification code**

The identification code includes the manufacturer, model code, serial number, and firmware revision levels and is sent in the following format: THORLABS,MMM,SSS,X.X.X

Where:       MMM is the model code  
              SSS is the serial number  
              X.X.X is the instrument firmware revision level

### **\*OPC – operation complete - set OPC bit**

### **\*OPC? – operation complete query – places a “1” in output queue**

When \*OPC is sent, the OPC bit in the Standard Event Register will set after all pending command operations are complete. When \*OPC? is sent, an ASCII “1” is placed in the Output Queue after all pending command operations are complete.

Typically, either one of these commands is sent after the INITiate command. The INITiate command is used to take the instrument out of idle in order to perform measurements. While operating within the trigger model layers, many sent commands will not execute. After all programmed operations are completed, the instrument returns to the idle state at which time all pending commands (including \*OPC and/or \*OPC?) are executed. After the last pending command is executed, the OPC bit and/or an ASCII “1” is placed in the Output Queue.

When \*OPC is sent, the OPC bit in the Standard Event Register will set after all pending command operations are complete. When \*OPC? is sent, an ASCII “1” is placed in the Output Queue after all pending command operations are complete.

### **\*RST – reset – return instrument to defaults**

When the \*RST command is sent, the instrument performs the following operations:

- Returns the instrument to the default conditions
- Cancels all pending commands.
- Cancels response to any previously received \*OPC and \*OPC? commands.

### **\*TST? – self-test query – run self test and read result**

Use this query command to perform the instrument self-test routine. The command places the coded result in the Output Queue. A returned value of zero (0) indicates that the test passed, other values indicate that the test failed.

### **\*WAI – wait-to-continue – wait until previous commands are completed**

The \*WAI command is a no operation command for the instrument and thus, does not need to be used. It is there for conformance to IEEE488.2.

## 2.2.3 PM200 specific SCPI Command Reference

See also SCPI Specification, Version 1999.0, May, 1999, <http://www.scpiconsortium.org> . All commands with a 'SCPI' checkmark are described in the SCPI specification.

All described commands work also with the PM100D, PM100A, PM100USB and PM160 instruments (with some limitations due to the hardware capabilities).

### 2.2.3.1 SYSTem subsystem commands

Command	Description																										
<b>SYSTem</b>	Path to SYSTem subsystem. (SCPI Vol.2 §21)																										
<b>:BEEPer</b>																											
<b>[ :IMMediate]</b>	Issue an audible signal. (SCPI Vol.2 §21.2.2)																										
<b>:STATe {ON 1 OFF 0}</b>	Activate/deactivate the beeper. (SCPI Vol.2 §21.2.3)																										
<b>:STATe?</b>	Return the state of the the beeper (SCPI Vol.2 §21.2.3)																										
<b>:ERRor</b>																											
<b>[ :NEXT]?</b>	Returns the latest <error code, "message">. (SCPI Vol.2 §21.8.8)																										
<b>:VERSion?</b>	Query level of SCPI standard (1999.0) . (SCPI Vol.2 §21.21)																										
<b>:DATE &lt;year&gt;,&lt;month&gt;,&lt;day&gt;</b>	Sets the instrument's calendar. (SCPI Vol.2 §21.7)																										
<b>:DATE?</b>	Query the instrument's calendar. (SCPI Vol.2 §21.7)																										
<b>:TIME &lt;hour&gt;,&lt;min&gt;,&lt;sec&gt;</b>	Sets the instrument's clock. (SCPI Vol.2 §21.19)																										
<b>:TIME?</b>	Query the instrument's clock. (SCPI Vol.2 §21.19)																										
<b>:LFRequency &lt;numeric value&gt;</b>	Sets the instrument's line frequency setting to 50 or 60Hz. (SCPI Vol.2 §21.13)																										
<b>:LFRequency?</b>	Query the instrument's line frequency setting. (SCPI Vol.2 §21.13)																										
<b>:SENSor</b>																											
<b>:IDN?</b>	<p>Query information about the connected sensor. This is a query only command. The response consists of the following fields: &lt;name&gt;, &lt;sn&gt;, &lt;cal_msg&gt;, &lt;type&gt;, &lt;subtype&gt;, &lt;flags&gt;</p> <table> <tr> <td>&lt;name&gt;</td><td>Sensor name in string response format</td></tr> <tr> <td>&lt;sn&gt;</td><td>Sensor serial number in string response format</td></tr> <tr> <td>&lt;cal_msg&gt;</td><td>calibration message in string response format</td></tr> <tr> <td>&lt;type&gt;</td><td>Sensor type in NR1 format</td></tr> <tr> <td>&lt;subtype&gt;</td><td>Sensor subtype in NR1 format</td></tr> <tr> <td>&lt;flags&gt;</td><td>Sensor flags as bitmap in NR1 format.</td></tr> </table> <table> <tr> <td><u>Flag:</u></td><td><u>Dec.value:</u></td></tr> <tr> <td>Is power sensor</td><td>1</td></tr> <tr> <td>Is energy sensor</td><td>2</td></tr> <tr> <td>Response settable</td><td>16</td></tr> <tr> <td>Wavelength settable</td><td>32</td></tr> <tr> <td>Tau settable</td><td>64</td></tr> <tr> <td>Has temperature sensor</td><td>256</td></tr> </table>	<name>	Sensor name in string response format	<sn>	Sensor serial number in string response format	<cal_msg>	calibration message in string response format	<type>	Sensor type in NR1 format	<subtype>	Sensor subtype in NR1 format	<flags>	Sensor flags as bitmap in NR1 format.	<u>Flag:</u>	<u>Dec.value:</u>	Is power sensor	1	Is energy sensor	2	Response settable	16	Wavelength settable	32	Tau settable	64	Has temperature sensor	256
<name>	Sensor name in string response format																										
<sn>	Sensor serial number in string response format																										
<cal_msg>	calibration message in string response format																										
<type>	Sensor type in NR1 format																										
<subtype>	Sensor subtype in NR1 format																										
<flags>	Sensor flags as bitmap in NR1 format.																										
<u>Flag:</u>	<u>Dec.value:</u>																										
Is power sensor	1																										
Is energy sensor	2																										
Response settable	16																										
Wavelength settable	32																										
Tau settable	64																										
Has temperature sensor	256																										

## 2.2.3.2 STATus subsystem commands

Command	Description
STATus	Path to STATus subsystem. (SCPI Vol.2 §20)
:MEASurement	Path to control measurement event registers
[:EVENTt]?	Read the event register
:CONDition?	Read the condition register
:PTRansition <value>	Program the positive transition filter
:PTRansition?	Read the positive transition filter
:NTRansition <value>	Program the negative transition filter
:NTRansition?	Read the negative transition filter
:ENABle <value>	Program the enable register
:ENABle?	Read the enable register
:AUXillary	Path to control measurement event registers
[:EVENTt]?	Read the event register
:CONDition?	Read the condition register
:PTRansition <value>	Program the positive transition filter
:PTRansition?	Read the positive transition filter
:NTRansition <value>	Program the negative transition filter
:NTRansition?	Read the negative transition filter
:ENABle <value>	Program the enable register
:ENABle?	Read the enable register
:OPERation	Path to control operation event registers
[:EVENTt]?	Read the event register
:CONDition?	Read the condition register
:PTRansition <value>	Program the positive transition filter
:PTRansition?	Read the positive transition filter
:NTRansition <value>	Program the negative transition filter
:NTRansition?	Read the negative transition filter
:ENABle <value>	Program the enable register
:ENABle?	Read the enable register
:QUEStionable	Path to control questionable event registers
[:EVENTt]?	Read the event register
:CONDition?	Read the condition register
:PTRansition <value>	Program the positive transition filter
:PTRansition?	Read the positive transition filter
:NTRansition <value>	Program the negative transition filter
:NTRansition?	Read the negative transition filter
:ENABle <value>	Program the enable register
:ENABle?	Read the enable register
:PRESet	Return status registers to default states.

### 2.2.3.3 DISPlay subsystem commands

Command	Description
<b>DISPlay</b>	Path to DISPlay subsystem. (SCPI Vol.2 §8)
<b>:BRIGhtness &lt;value&gt;</b>	Set the display birghtness. (SCPI Vol.2 §8.2)
<b>:BRIGhtness?</b>	Return the display birghtness value. (SCPI Vol.2 §8.2)
<b>:CONTRast &lt;value&gt;</b>	Set the display contrast. (SCPI Vol.2 §8.4)
<b>:CONTRast?</b>	Return the display conrast value. (SCPI Vol.2 §8.4)

### 2.2.3.4 CALibration subsystem commands

Command	Description
<b>CALibration</b>	Path to CALibration subsystem. (SCPI Vol.2 §5)
<b>:STRing?</b>	Returns a human readable calibration string. This is a query only command. The response is formatted as string response data.

### 2.2.3.5 SENSE subsystem commands

Command	Description
<b>SENSe</b>	Path to SENSE subsystem. (SCPI Vol.2 §18)
<b>AVERage</b>	
<b>[ :COUNT ] &lt;value&gt;</b>	Sets the averaging rate (1 sample takes approx. 3ms)
<b>[ :COUNT ] ?</b>	Queries the averaging rate
<b>CORRection</b>	
<b>[ :LOSS [ :INPut [ :MAGNitude ] ] ] { MINimum   MAXimum   DEFault   &lt;numeric_value&gt; }</b>	Sets a user attenuation factor in dB
<b>[ :LOSS [ :INPut [ :MAGNitude ] ] ] ? { MINimum   MAXimum   DEFault }</b>	Queries the user attenuation factor
<b>COLLect</b>	
<b>ZERO</b>	
<b>[ :INITiate ]</b>	Performs zero adjustment routine
<b>ABORt</b>	Aborts zero adjustment routine
<b>STATe?</b>	Queries the zero adjustment routine state
<b>MAGNitude?</b>	Queries the zero value
<b>BEAMdiameter { MINimum   MAXimum   DEFault   &lt;numeric_value&gt; [mm] }</b>	Sets the beam diameter in mm
<b>BEAMdiameter? [ { MINimum   MAXimum   DEFault } ]</b>	Queries the beam diameter
<b>WAVelength { MINimum   MAXimum   &lt;numeric_value&gt; [nm] }</b>	Sets the operation wavelength in nm
<b>WAVelength? [ { MINimum   MAXimum } ]</b>	Queries the operation wavelength
<b>POWer</b>	
<b>[ :PDIODE ]</b>	Sets the photodiode response value in A/W

Command	Description
<code>[ :RESPonse] {MINimum MAXimum DEFAULT &lt;numeric_value&gt;[A]}</code>	
<code>[ :RESPonse]? [{MINimum MAXimum DEFAULT}]</code>	Queries the photodiode response value
<code>:THERmopile</code>	
<code>[ :RESPonse] {MINimum MAXimum DEFAULT &lt;numeric_value&gt;[V]}</code>	Sets the thermopile response value in V/W
<code>[ :RESPonse]? [{MINimum MAXimum DEFAULT}]</code>	Queries the thermopile response value
<code>ENERgy</code>	
<code>[ :PYRO]</code>	
<code>[ :RESPonse] {MINimum MAXimum DEFAULT &lt;numeric_value&gt;[V]}</code>	Sets the pyro-detector response value in V/J
<code>[ :RESPonse]? [{MINimum MAXimum DEFAULT}]</code>	Queries the pyro-detector response value
<code>CURRent [:DC]</code>	
<code>RANGE</code>	
<code>AUTO {OFF 0 ON 1}</code>	Switches the auto-ranging function on and off
<code>AUTO?</code>	Queries the auto-ranging function state
<code>[ :UPPer] {MINimum MAXimum &lt;numeric_value&gt;[A]}</code>	Sets the current range in A
<code>[ :UPPer]? [{MINimum MAXimum}]</code>	Queries the current range
<code>REFerence {MINimum MAXimum DEFAULT &lt;numeric_value&gt;[A]}</code>	Sets a delta reference value in A
<code>REFerence? [{MINimum MAXimum DEFAULT}]</code>	Queries the delta reference value
<code>STATe {OFF 0 ON 1}</code>	Switches to delta mode
<code>STATe?</code>	Queries the delta mode state
<code>ENERgy</code>	
<code>RANGE</code>	
<code>[ :UPPer] {MINimum MAXimum &lt;numeric_value&gt;[J]}</code>	Sets the energy range in J
<code>[ :UPPer]? [{MINimum MAXimum}]</code>	Queries the energy range
<code>REFerence {MINimum MAXimum DEFAULT &lt;numeric_value&gt;[J]}</code>	Sets a delta reference value in J
<code>REFerence? [{MINimum MAXimum DEFAULT}]</code>	Queries the delta reference value
<code>STATe {OFF 0 ON 1}</code>	Switches to delta mode
<code>STATe?</code>	Queries the delta mode state
<code>FREQuency</code>	
<code>Range</code>	

Command	Description
[UPPer]?	Queries the frequency range
LOWer?	
POWer[:DC]	
RANGe	
AUTO {OFF 0 ON 1}	Switches the auto-ranging function on and off
AUTO?	Queries the auto-ranging function state
[:UPPer] {MINmum MAXimum <numeric_value>[W]}	Sets the current range in W
[:UPPer]? [{MINimum MAXimum}]	Queries the current range
REfERENCE {MINimum MAXimum DEfault <numeric_value>[W]}	Sets a delta reference value in W
REfERENCE? [{MINimum MAXimum DEfault}]	Queries the delta reference value
STATe {OFF 0 ON 1}	Switches to delta mode
STATe?	Queries the delta mode state
UNIT {W DBM}	Sets the power unit W or dBm
UNIT?	Queries the power unit
VOLTage[:DC]	
RANGe	
AUTO {OFF 0 ON 1}	Switches the auto-ranging function on and off
AUTO?	Queries the auto-ranging function state
[:UPPer] {MINmum MAXimum <numeric_value>[V]}	Sets the current range in V
[:UPPer]? [{MINimum MAXimum}]	Queries the current range
REfERENCE {MINimum MAXimum DEfault <numeric_value>[V]}	Sets a delta reference value in V
REfERENCE? [{MINimum MAXimum DEfault}]	Queries the delta reference value
STATe {OFF 0 ON 1}	Switches to delta mode
STATe?	Queries the delta mode state
PEAKdetector	
[:THReshold] {MINimum MAXimum DEfault <numeric_value>}	Sets the trigger level in % for the energy mode
[:THReshold]? [{MINimum MAXimum DEfault}]	Queries the trigger level setting

### 2.2.3.6 INPut subsystem commands

Command	Description
<b>INPut</b> [:PDIode] :FILTER [:LPASs] [STATE] {OFF 0 ON 1} [STATE]?	Sets the bandwidth of the photodiode input stage Queries the bandwidth of the photodiode input stage
:THERmopile :ACcelerator [STATE] {OFF 0 ON 1} [STATE]? :AUTO {OFF 0 ON 1} ? :TAU {MINimum  MAXimum DEFAULT  <numeric_value>[s]} :TAU? [{MINimum  MAXimum DEFAULT}]	Sets the thermopile accelerator state Queries the thermopile accelerator state Sets the thermopile accelerator to auto mode Queries thermopile accelerator auto mode Sets thermopile time constant 0-63% in s  Queries the thermopile time constant in s
:ADAPter [:TYPE] {PHOTodiode  THERmal PYRo} [:TYPE]?	Sets default sensor adapter type  Queries default sensor adapter type

### 2.2.3.7 Measurement commands

Command	Description
INITiate[:IMMediate]	Start measurement
ABORt	Abort measurement
CONFigure[:SCALar]	
[:POWer]	Configure for power measurement
:CURRent[:DC]	Configure for current measurement
:VOLTag[:DC]	Configure for voltage measurement
:ENERgy	Configure for energy measurement
:FREQuency	Configure for frequency measurement
:PDENsity	Configure for power density measurement
:EDENsity	Configure for energy density measurement
:RESistance	Configure for sensor presence resistance measurement
:TEMPerature	Configure for sensor temperature measurement
MEASure[:SCALar]	
[:POWer]	Performs a power measurement
:CURRent[:DC]	Performs a current measurement
:VOLTag[:DC]	Performs a voltage measurement
:ENERgy	Performs a energy measurement
:FREQuency	Performs a frequency measurement
:PDENsity	Performs a power density measurement
:EDENsity	Performs a energy density measurement
:RESistance	Performs a sensor presence resistance measurement
:TEMPerature	Performs a sensor temperature measurement
FETCh?	Read last measurement data (SCPI Vol.2 §3.2)
READ?	Start new measurement and read data (SCPI Vol.2 §3.3)
CONFigure?	Query the current measurement configuration.